

P. Merkl^a, T. Braatz^a, M. J. Prandolini^{a,b,*}, H. Goudarzi^a, S. Starosielec^a, T. Golz^a, B. Manschwetus^a, M. Schulz^a, and R. Riedel^a

^aClass 5 Photonics GmbH, Notkestr. 85, 22607 Hamburg, Germany
^bUniversität Hamburg, Institut für Experimentalphysik, Luruper Chaussee 149, 22761 Hamburg, Germany
*mark.prandolini@class5photonics.com


Motivation

Class 5 Photonics delivers ultrafast lasers for state-of-the-art applications, such as pump-probe-spectroscopy [1]/ -microscopy [2], attosecond science [3], particle acceleration [4], advanced micromachining and fabrication [5]. These laser systems require both high average power and high peak power. Many applications, in addition to the high average and peak power, require wavelength tunability, which can be provided by optical parametric chirped-pulse amplifier (OPCPA) technology. With a focus on customer demands, Class 5 Photonics has developed a tunable, high power OPCPA with sub-50 fs pulses ranging from the ultraviolet to the near-infrared spectral regions, which we call the *White Dwarf HE (OPCPA)*.


An important advantage of our OPCPA design is that it can be pump by various commercially available Yb-doped pump lasers at a fixed wavelength of 1030 nm, including the Amplitude *Tangor*, Trumpf Scientific *Dira* and the Light Conversion *Carbide*.

In this poster, we compare the system performances in terms of conversion efficiency, beam and pulse quality and show how the pump parameters translate into OPCPA output. Further details can be found in reference [6].

System Layout




Performance: Tunability from UV to NIR




Module	Parameter	Measurements White Dwarf HE OPCPA		
Pump laser	Model	Light conversion Carbide CB3-80W	Amplitude Tangor 300	Trumpf Scientific DIRA 500-10
	Pump power	20 W	150 W	30 W
	Repetition rate	10 kHz	100 kHz	10 kHz
	Pulse energy	2 mJ	1.2 mJ	3 mJ
	Pulse duration	350 fs	500 fs	800 fs
NIR OPCPA	Spectral tuning range	700 – 900 nm		
	Output power	> 1.5 W	> 12 W	> 2 W
	Pulse energy	> 150 μJ	> 120 μJ	> 230 μJ
	Power stability	< 0.5 % rms	< 0.3 % rms	< 2 % rms
	Pulse duration	< 30 fs	< 35 fs	< 35 fs
SH Module	M2 (x,y)	< 1.4 in both axes	< 1.5 in both axes	< 1.1 in both axes
	Spectral tuning range	350 – 450 nm		
	Output power	> 0.4 W	> 1.5 W	> 0.6 W
	Pulse energy	> 40 μJ	> 15 μJ	> 60 μJ
	Power stability	-	-	< 3 % rms
TH Module	Pulse duration	< 50 fs	< 50 fs	< 60 fs
	Spectral tuning range	235 – 300 nm		
	Output power	> 90 mW	> 200 mW	> 150 mW
	Pulse energy	> 9 μJ	> 2 μJ	> 15 μJ
	Power stability	-	-	< 3 % rms
FH Module	Pulse duration	< 100 fs (without prism comp.)	< 40 fs (with prism comp.)	< 100 fs (without prism comp.)
	Spectral tuning range	-	200-220 nm	-
	Output power	-	1 - 10 mW	-
	Pulse energy	-	10 - 100 nJ	-

Summary of laser output parameters, pumped by three different commercially available Yb-pump lasers

Autocorrelation: (a) 800, (b) 400 and (c) 266 nm



Beam Stability and Profile



References

[1] F. Mahmood, et al, *Nature Phys.* **12**, 306-310 (2016); J. Oh, et al, *J. Phys. Chem C* **127**, 5004-5012 (2023)
[2] S. Weisenburger, et al, *Cell* **177**, 1050 – 1066 (2019); A. Klioutchnikov, et al, *Nature Methods*, **20**, 610-616 (2023)
[3] D. Hui, et al. *Sci. Adv.* **10**, eadp5805 (2024)
[4] T. Chlouba, et al, *Nature* **622**, 476 – 480 (2023)
[5] C. Porter, et al, "Soft x-ray: novel metrology for 3D profilometry and device pitch overlay", *Proc. SPIE 12496, Metrology, Inspection, and Process Control XXXVII* (2023)
[6] V. Shumakova, et al, "Tunable, high-power deep UV to NIR source of femtosecond pulses utilizing various Yb-doped pump-laser architectures" *Proc. SPIE 13347 133470J-1*

